Morriston Hospital, Morriston, Swansea, SA6 6NL Tel: (01792) 702222 Fax: (01792) 703632

HTN: 01789 100

Ysbyty Treforys, Treforys, Abertawe, SA6 6NL Ffôn: (01792) 702222 Ffacs: (01792) 703632

HTN: 01789 100





## Burns and Plastic Surgery Occupational Therapy

EB 838

## Instructions for using Silicone Sheeting

Your Product is
Tour Frounct is

## **Directions for use**

As instructed by your key worker – initially for 2-3 hours per day, then over a period of 1-2 weeks gradually increase by 1 or 2 hours daily, allowing for washing periods.

- Open the peel pack and remove the sheet
- Cut (if necessary) to appropriate size overlapping the scar by a minimum of one centimetre.
- When applying the sheet, make sure the area is clean and dry and free from cream (remove any excess cream using a soft tissue).
- Remove the protective film and apply the silicone sheeting over the scar. Avoid stretching when applying over a joint.
- Remove the sheet three times a day for inspection, washing of the skin and scar
  massage unless otherwise instructed by your key worker. The silicone sheet can then
  be re-applied once the area is clean and dry and free from cream.

## Care of sheeting

- When not in use, the silicone sheeting should be returned to its plastic backing.
- The sheet should be stored in a cool place away from direct heat and sunlight.
- Wash gel daily in lukewarm water and mild soap. Air dry and re-apply.
- When adherent properties of the sheet are no longer sufficient discard and replace with new piece. Each piece should last a minimum of 2 weeks (1 week with Mepiform sheet) and up to a maximum of 4 weeks.

If you require any further guidance or have any questions about the silicone sheeting or if the scar area shows any evidence of skin irritation, rash, open or weeping areas, remove gel and contact your key worker.

Telephone: (01792) 703980 (Answer machine available)